



Container Closure Integrity Testing  
for the  
Pharmaceutical Industry

Tobias Schappeler  
Managing Director

# About us

---



- Established in 1989, family owned
- Sound experience in chemistry & pharma
- 21 Team Members, growing quickly
  - Pharmaceutical Chemists (some with validation experience)
  - Analytical Chemists
  - Engineers
  - Equipment Specialists
- Covering Australia and New Zealand with offices in Sydney, Brisbane, Melbourne, Perth, NZ with Partner
- Services include equipment consulting, application and technology support, equipment sales, installation, training and after sales service including recalibration services


# What is it?

---

## Container Closure Integrity Testing

- (US FDA) obligation for “pharma” is that they must ensure *“the container-closure system to maintain the integrity of its microbial barrier, and hence the sterility of a drug product throughout its shelf life”*
  
- The key **sources** for ingress contamination are
  - **humidity (H<sub>2</sub>O)**
  - **oxygen and**
  - **Microbiological**
  
- Container integrity is also important to retain volatiles in your product
  
- Key guideline describing different test methods and selection process for sterile products is **<USP 1207>** published by the PDA

# Microbial Ingress?



- Based on Lee Kirsch studies (presented in PDA J 51.5, 1997 p 195-202) the **critical leak size for microbiological ingress is 0.2  $\mu\text{m}$**
- The risk of contamination increases with defect size
- With a defect of 0.7 $\mu\text{m}$  the risk of contamination is already >60%
- Above 5 $\mu\text{m}$  the risk of contamination is >80%
- Only 2 test current methods available can detect <0.2  $\mu\text{m}$

# Probabilistic Testing

## Blue Dye Ingress Testing

- the chance to detect a 10 $\mu$ m leak is only about 70%
- defects below 5 $\mu$ m are most likely not detected
- At 5 $\mu$ m our risk of contamination is >80%
- 0.2  $\mu$ m is the ultimate target we aim for

The detection limit of Blue Dye Tests varies with:


Technical:

Leak size, type, length, material of construction, type of blockage, tracer concentration, surface tension, cleanliness, tracer compatibility with product or immersion fluid

Operational:

Ambient temperature, ambient pressure, sample positioning, inspection conditions, operator training/skill, sample preparation

## USP/PhEur Dye Ingress Test Samples



## Common “Technology Based” Test Methods

Leak Test Method	Measurement Outcome	Detection range
<b>Tracer-gas</b> (Helium Mass Spectrometry)	Helium flow (mbar.l/s)	<0.1µm to 10µm
<b>Laser-Headspace</b> (Frequency modulated spectroscopy)	[O <sub>2</sub> ] and/or [CO <sub>2</sub> ] (%) Gas pressure	<0.1µm to >50µm
<b>AMI*</b> (Optical Emission Spectroscopy)	Leakage (N <sub>2</sub> , Ar, CO <sub>2</sub> , H <sub>2</sub> O..) (mbar.l/s)	<0.5µm to >50µm
<b>Mass Extraction</b> (Micro/Mass flow sensors)	Mass Flow (µg/min)	>1.0µm to >50µm
<b>HVLD</b> (Leakage current)	Electrical current (µA)	>1.0µm to >50µm
<b>Pressure Decay</b>	Pressure drop (mbar/s)	>1.0µm to >50µm
<b>Vacuum Decay</b>	Pressure rise (mbar/s)	>1.0µm to >50µm

\*not yet recognized in USP 1207, but has been presented on PDA conferences – emerging technology

# 3 Leak Detection Solutions

## MICRO-FLOW AND MASS EXTRACTION

Air micro-flow sensor



ASTM F3287-17

## HELIUM MASS SPECTROMETRY

Magnetic deflection spectrometer



ASTM F2391-05




## OPTICAL EMISSION SPECTROMETRY

Multi-gas analyser (N<sub>2</sub>, CO<sub>2</sub>, Ar, H<sub>2</sub>O)




3 different technologies for CCIT solutions

→ because there is no one size fits all solution

-  Non destructive test options
-  Applicable for non-porous containers
-  Global pass / fail test

# Helium Leak Detection

- High Sensitivity and Quantitative
  - Mass Spectrometer (magnetic deflection)
  - Down to  $10^{-10}$  mbar.l/s (sub-micron orifice)
- High Selectivity
  - Low natural background of 5 ppm (in air)
    - Sensitivity can be affected by background accumulation
  - High permeability, diffusivity & solubility
- Helium flows through cracks
  - Much smaller and faster than air





# ASM 2000

- Turn-key equipment dedicated for pharma
  - Based on high performance helium leak detector
  - All in one, including helium charging module
  - PLC and HMI (3",5 touch screen)
  - Customized test tooling according to the part to test
  - Trolley includes all vacuum pumps
  - Data storage / 3 access level / PDF test reports
  - CFR 21 P 11 soon to be released
- *MALL TEST = Maximum allowable Leak Limit*
- *Global Pass & Fail*



# R&D Testing with ASM 2000

## Experimental Glass Bottle Cap Test - how to proceed ?

- Helium injection inside the bottle must be controlled & performed during the test sequence
- Leak testing must be performed before helium permeation appears (plastic)

### ASM 2000 Test sequence:

- 1/. Start test on ASM 2000
- 2/. Bottle evacuation (remove air)
- 3/. Helium Charging (Patm)
- 4/. Helium test
- 5/. Helium evacuation
- 6/. N<sub>2</sub> venting/purge
- 7/. Stop test on HLD / Venting

Fully automated process, started with one button, settings can be customized

A 6mm diameter hole has been drilled in the bottom of the bottles




## In production = non destructive test approaches

-  Injection of He tracer gas prior to container closure
-  “Bombing” Test – apply pressurized He

# Non Destructive Helium Testing


**Container sealed**



1) Bombard with helium tracer gas whilst sealed

2) Global vacuum pass & fail (MALL) test

**Container open / unsealed for filling**




1) Inject helium tracer gas before sealing the package

2) Global vacuum pass & fail (MALL) test

# Optical Emission Spectroscopy

A Plasma is the 4<sup>th</sup> fundamental state of matter = hot ionized gas (instable).


Plasma de-excitation gives light emission.



Test chamber ( $P < 5 \cdot 10^{-2}$  mbar)


Sealed component  
( $N_2$ , Ar,  $CO_2$ ...)

$Q_{N_2}$   
( $atm \cdot cm^3/s$ )



200nm

850nm



patented OES multi-gas "ALPS" sensor

# AMI 1000

## Detection limits / cycle times for different package types

Samples	Sensitivity Orifice diameter <sup>3)</sup> Air/N <sub>2</sub> Leakage		Test duration	Advantages
	Air/N <sub>2</sub> leak	Water leak		
Blisters	0.4 µm 2 · 10 <sup>-5</sup> mbar l/s	n.a.	> 20-30 sec	Highest sensitivity test method available on the market Outgazing of the drug itself can be used for gross leak detection Applicable to peeling blisters
Syringes & Vials	0.4 µm 2 · 10 <sup>-5</sup> mbar l/s	2 µm	> 15 sec	Air and water detected simultaneously Test per batch to increase the throughput
	0.2 µm 6 · 10 <sup>-6</sup> mbar l/s		~45 sec. (high sensitivity mode)	MALL level can be achieved in high sensitivity mode
IV bags	0.4 µm 2 · 10 <sup>-5</sup> mbar l/s	3 µm	> 20 sec	Air and water detected simultaneously
Plastic bottles	0.5 µm 4 · 10 <sup>-5</sup> mbar l/s	n.a.	> 20 sec	Test per batch (up to 50 or 100) to increase the throughput


<sup>3)</sup>Sharp edge orifice as defined in USP <1207> guidelines

- OES is unique amongst leak detection due to its ability to track various gases, moisture etc (N<sub>2</sub>, H<sub>2</sub>O, Ar, CO<sub>2</sub>...)
- Detection limits depends on tracked gas species
- All can be tracked simultaneously


# Mass Extraction

- Measurement of mass flow rate from test chamber to a reference
- Flow equals defect (down to 1  $\mu\text{m}$ )
- Gas based – typically air or nitrogen
- For vacuum testing water vapour flow is used


Mass Conservation law:  
Mass extracted =  
mass leaked at  
steady state




**1 Standby** – Vacuum Reservoir:  $P_0$   
– IGLS: no flow



**2 Fill** – All branches:  $P_0$   
– IGLS: no flow  
– UUT: leaks



**3 Stabilize** – Leak: increases:  $P_{\text{chamber}}$   
– IGLS: begins flow



**4 Test** – Steady flow thru leak  
– IGLS: measure flow thru leak


# Example IV Bag Testing

Storage area for  
laptop,  
packaging



ATC  
calibrated  
leak device  
(inside)

Vacuum  
System



Fillers for  
different bag  
sizes

Customized  
test chamber  
(IV-bag)

ME2 Integrity  
Test Unit


# Example IV Bag Testing

Storage area for  
laptop,  
packaging



ATC  
calibrated  
leak device  
(inside)

Vacuum  
System



Fillers for  
different bag  
sizes

Customized  
test chamber  
(IV-bag)

ME2 Integrity  
Test Unit




# ASTM F3287 – 17 (Mass Extraction) Result Extract



- Glass vials and LDPE Bottles Mass Extraction tests detected 1µm and 2µm defects at all labs and samples at over 95% confidence level
- Glass syringes Mass Extraction tests detected 1µm air filled syringes and 2µm air and water filled syringes at all labs and samples. 2µm were detected at a confidence level equal or greater that 95%
- 1µm liquid filled syringe plugged – suspected by silicon lubricant

	Package Description	Sample Qty.	Qty. of Tests	Qty. of Failed Tests	Qty. of Passed Tests	Success %
Glass Vial 2 ml	Liquid Filled – Negative Control	10	120	0	120	100 %
	Air Filled – Negative Control	10	120	0	120	100 %
	1 µm micropipette – Liquid Filled	3	36	36	0	100 %
	1 µm micropipette – Air Filled	3	36	36	0	100 %
	2 µm micropipette – Liquid Filled	3	36	36	0	100 %
	2 µm micropipette – Air Filled	3	36	36	0	100 %
	5 µm micropipette – Liquid Filled	3	36	36	0	100 %
	5 µm micropipette – Air Filled	3	36	36	0	100 %
	10 µm micropipette – Air Filled	3	36	36	0	100 %
	LDPE Bottle 4 ml	Liquid Filled – Negative Control	10	120	0	120
Air Filled – Negative Control		10	120	0	120	100 %
1 µm micropipette – Liquid Filled		3	36	36	0	100 %
1 µm micropipette – Air Filled		3	36	36	0	100 %
2 µm micropipette – Liquid Filled		3	36	36	0	100 %
2 µm micropipette – Air Filled		3	36	36	0	100 %
5 µm micropipette – Liquid Filled		3	36	36	0	100 %
5 µm micropipette – Air Filled		3	36	36	0	100 %
10 µm micropipette – Air Filled		3	36	36	0	100 %
Glass Syringe 1 ml		Air Filled – Negative Control	10	120	0	120
	1 µm micropipette – Air Filled	3	36	36	0	100 %
	2 µm micropipette – Air Filled	3	36	36	0	100 %
	5 µm micropipette – Air Filled	3	36	36	0	100 %
	10 µm micropipette – Air Filled	3	36	36	0	100 %
	Liquid Filled – Negative Control	10	120	0	120	100 %
	1 µm micropipette – Liquid Filled	3	36	0	36	0 %
	2 µm micropipette – Liquid Filled	3	36	36	0	100 %
5 µm micropipette – Liquid Filled	3	36	36	0	100 %	

# Example IV Bags- Inline Testing



# Example Vial Batch



# Summary

- Scitek offers highly reliable DETERMINISTIC test methods for different kinds of pharmaceutical packages and drug types (liquid or solid) – applicable for...



- We offer project support including design of customized tooling, IQ/OQ support as well as FAT / SAT support
- All our tools are factory calibrated based on traceable leak standards
- CFR 21 part 11 compliant software is available for He and Mass Extraction testing

# Whitepaper

## Overview of existing CCIT Technologies (not all covered today)

	Helium Mass Spectrometry	O.E.S (Optical Emission Spectroscopy)	Mass Extraction	Vacuum Decay	Viscon (Deflexion)	HSA (Head Space Analysis)	HVLD (High Voltage)	Dye Ingress	Microbial Challenge	
Deterministic	Yes (Yes)	Yes	Yes	Yes	Yes	Yes	Yes			
Non-Destructive	only for open containers	Yes	Yes	Yes	Yes	Yes	Yes			
Quantitative	Yes	Yes	Yes	Yes		Yes				
Sample preparation	He charging Plausibility test.	No sample preparation				Storage time	No sample preparation	Immersion in dye or microbial media		
Test pressure		Vacuum			Atmospheric Pressure		Shallow Vacuum	Atmospheric Pressure		
Detection range (Sharp edge orifice)	0,01 < Q < 10 µm	> 0,2 µm	> 1 µm	> 5 µm	> 5 µm	> 0,01 µm	10-40 µm	> 20 µm	> 0,2 µm	
Drug Product Limitations		Lyophilized (dry) or liquid drugs Plugging risk for small defects for protein based drugs Container must handle 1 bar differential pressure				Lyophilized drugs	Conductive liquid drugs	Light colored drugs	-	
Container Limitations	He Permeation	High outgassing			Container Design (Semi-rigid or flexible)	Rigid & Transparent	Non-conductive material	Non-porous material		
Method Limitations	Require gas headspace or liquid inside the container									
	<ul style="list-style-type: none"> <li>Difficult to set-up</li> <li>Requires proper He gas management</li> <li>Requires plausibility test to valid the test result.</li> </ul> <p>Not practical for mass production testing</p>	Outgassing of the container and the drug type will impact the test duration and the detection limit			Sensitivity depends on the product design: <ul style="list-style-type: none"> <li>Headspace volume</li> <li>Size of the cavity</li> <li>Shape of the container</li> </ul>		Requires waiting time before actual testing (hours up to weeks)	Test only at the point of electrode contact, with liquid behind	Destructive	Long (few weeks) and Expensive
		Detection limit is depending on packaging and drug type	Free volume inside the test chamber can limit sensitivity -> Test chamber must be optimized for each format parts.		Requires positive control to calibrate the equipment		Waiting time depends on the gas headspace and detection limit.	Limited usage for flexible packaging's.	Probabilistic	
	Detection limit depends on the gas used for the detection	Sensitive to temperature and/or volume variations				Headspace needs to be either vacuum or 100% Nitrogen	No real quantitative measurement	Poor sensitivity		
Method Advantages	<ul style="list-style-type: none"> <li>High selectivity (He)</li> <li>High sensitivity test</li> </ul> <p>Possibility to localize the leak position with sniffing.</p>	<ul style="list-style-type: none"> <li>Selectivity: can detect simultaneously gas species (N<sub>2</sub>, H<sub>2</sub>O, Ar, CO<sub>2</sub>,...)</li> </ul> <p>Can test multiple containers with high sensitivity at the same time.</p>	<ul style="list-style-type: none"> <li>High sensitivity detection of water leakage</li> <li>Robust technology</li> </ul>	Simple	<ul style="list-style-type: none"> <li>Identification of the leaky cavity or container.</li> </ul> <p>Can test multiple containers with high sensitivity at the same time.</p>	<ul style="list-style-type: none"> <li>High selectivity (O<sub>2</sub>)</li> <li>Very fast, high throughput can be achieved</li> </ul>	Very fast, high throughput can be achieved	Low cost equipment	Direct measurement of the biological contamination	
Comments	Mainly used for the design and qualification phase of the packaging's, not practical for mass production testing.	Highly versatile and sensitive test for different drug / packaging systems Can be used in laboratory or as IPC in production.	Highly versatile and sensitive test for different drug / packaging systems In-line option available.	Older production test method. Reduced reliability for measurements at limit of detection.	Mainly used for blister packs.	Indirect leak test, we measure the consequence of oxygen ingress through defects.	Very fast method for production test, limited usage for flexible packaging's.	Easy to understand		Widely-used for decades Industry & regulatory familiarity

# Thank you

---



## Scitek Contacts

- [scitek.com.au](http://scitek.com.au)
- [medicinalcannabisproduction.com.au](http://medicinalcannabisproduction.com.au)
- 1800 023 467
- [sales@scitek.com.au](mailto:sales@scitek.com.au)

## Personal Contacts

- [tobias@scitek.com.au](mailto:tobias@scitek.com.au)
- 0413 756 224